

Paper Pulp Bioplastic

(This is a great recipe with very strong material characteristics!!)

Material Properties:

- Soft leather feel • Flexible and elastic
- Translucent • Easy to manipulate, cut etc.
- 500 ml water
- 11g Agar powder
- 19g Vegetable glycerin
- 8g cardboard pulp (or more!!!)

Bring water to a boil, add glycerine and stir until dissolved, add agar powder stirring slowly for 40 minutes, cut recycled paper into a blender and add water then blend, add agar mixture into paper pulp mixture and blend together, leave to dry on a flat surface.